
Grapes

Several Pounds of grapes are needed to fill 4 medium trays

Adventures in Freeze Drying for Beginners: The Grape Experiment

Ingredients:





Directions:

- 1. Use your favorite method for cleaning grapes, a chemical free way is to dip in boiling water and then rub each grape with a paper towel.. You can also add a cap full of vinegar to a sink of cool water and let the grapes soak, make sure to rinse thoroughly.
- 2. Slice each grape in half
- 3. Line trays with parchment or silicone if you choose, easier to remove from trays if they are lined
- 4. Spread out in a layer
- 5. Freeze

Cycle Time: My cycle time was ~24 hours in a a Medium Freeze Dryer with a Premier Pump	Rehydration: Not intended for rehydration
---	---

